

From venous blood to VAMS – adaptation of a commercial CE-IVD assay kit for the determination of Tacrolimus, Sirolimus, Everolimus and Cyclosporin A.

Rafal Szewczyk^{1,4}; Arkadiusz Kocur^{2,3}; Agnieszka Czajkowska³; Karina Izdebska⁴; Maciej Stopa^{1,4}. 1 - LabExperts sp. z o.o., Gdańsk, Poland; 2 - Department of Drug Chemistry, Pharmaceutical and Biomedical Analysis, Faculty of Pharmacy, Medical University of Warsaw, Warsaw, Poland; 3 - Therapeutic Drug Monitoring, Clinical Pharmacokinetics and Toxicology Laboratory Unit, Department of Clinical Biochemistry, The Children's Memorial Health Institute, Warsaw, Poland; 4 - Bioanalytic Sp. z o.o., Gdańsk, Poland

BACKGROUND

TDM of immunosuppressants is extremely important topic, especially when it comes to sample collection at home coupled with simple and fast sample preparation and analysis. The study is focused on the check if optimization of ready-to-use kit (Chromsystems CE-IVD kit: 93000) for determination of Tacrolimus (TAC), Sirolimus (SIR), Everolimus (EVR) and Cyclosporin A (CSA) in venous blood can be successfully applied for VAMS-collected blood and preliminary venous blood to Mitra® (VB/M) correlation study.

METHODS

Venous blood (Chromsystems procedure)

Mitra® - optimized procedure based on Chromsystems kit



- Pipette 50 µL of sample/calibrator/Mass-Check® into a 1.5 mL tube.
- Add 25 µL of the reconstituted ISTD solution.
- Add 100 µL of extraction buffer, vortex for 10 s and incubate for 2 minutes at 20-25 °C.
- Add 250 µL of precipitating reagent, vortex for 1 minute and incubate for 2 minutes at 20-25 °C.
- Centrifuge for 5 minutes at 15000 x g.
- Transfer the supernatant to the vial.
- Inject 2-5 µL (SCIEX 5500+ QTRAP, Exion AC LC-MS/MS).

- Pipette min. 40 µL sample/calibrator/Mass-Check® on petri dish.
- Collect the material on the Mitra® sampler (10 µL).
- Dry samples at least 2h in the dark, 20-25°C.
- Transfer the Mitra® sampler to a 2 mL tube.
- Prepare an extraction buffer containing ISTD (20 µL ISTD + 100 µL extraction buffer x number of samples + 1)
- Add 120 µL of extraction buffer containing ISTD.
- Vortex for 10 s and then sonicate for 5 min at 20-25 °C.
- Add 80 µL of precipitating reagent, vortex for 10 s and sonicate for 5 min at 20-25 °C.
- Centrifuge for 5 minutes at 15000 x g.
- Transfer 100 µL of supernatant to the insert (100 µL on the feet) placed in the vial.
- Inject 10 µL (SCIEX 5500+ QTRAP, Exion AC LC-MS/MS).



Q1	Q3	Dwell time [ms]	Id	DP	EP	CE	CXP	IS Parameters
821.5	768.5	50	Tacrolimus	110	10	30	12	CUR 30
824.5	771.5	50	Tacrolimus 13CD2	110	10	30	12	CAD 9
931.6	864.5	50	Sirolimus	110	10	25	14	IS 5500
934.6	864.5	50	Sirolimus-D3	110	10	25	14	TEM 300
975.6	908.6	50	Everolimus	110	10	27	14	GS1 45
979.6	912.6	50	Everolimus-D4	110	10	27	14	GS2 45
1219.9	1202.8	50	Cyclosporin A	110	10	26	16	
1223.9	1206.8	50	Cyclosporin A-D4	110	10	26	16	

Data processing, quantitation and basic statistical analysis was done in SciexOS 3.3 software (SCIEX).

RESULTS

Tab. 1. Validation results for collected on 10 µL Mitra® sampler lowest analytes concentrations provided with the kit (CAL 1, 4 replicates). LOD was not determined, however additional analyses showed that lower than 50-100 µg/ml concentration for each compound can be achieved.

Compound	LLOQ				ULOQ [ng/ml]	Linearity [R ≥ 0.995]
	C [ng/ml]	S/N ≥ 10	Accuracy [100±15%]	%CV [±15%]		
TAC	1.28	11.66	102.75	12.1	37.50	0.99907
SIR	2.10	15.15	100.05	12.3	43.40	0.99872
EVR	2.20	14.82	100.95	8.7	42.25	0.99826
CSA	24.20	52.68	111.30	7.9	772.00	0.99928

Signal-to-noise (S/N) was determined with peak-to-peak algorithm and averaged from 4 replicates.

Tab. 2. Results of preliminary study for VB/M correlation included 15 patients per compound where venous blood and 10 µL Mitra® samples were collected at the same time. No patients for EVR were tested during the study.

TAC				SIR				CSA			
Sample	Venous blood [ng/ml]	Mitra [ng/ml]	VB/M	Sample	Venous blood [ng/ml]	Mitra [ng/ml]	VB/M	Sample	Venous blood [ng/ml]	Mitra [ng/ml]	VB/M
1	4.011	5.257	0.763	157	2.196	2.771	0.792	187	56.379	44.22	1.275
2	9.085	7.686	1.182	158	3.798	3.026	1.255	188	92.872	78.181	1.188
3	6.65	12.425	0.535	159	1.089	1.867	0.583	189	81.334	71.476	1.138
4	5.198	8.861	0.587	160	3.763	5.619	0.670	190	175.848	174.666	1.007
5	7.579	9.945	0.762	161	2.212	3.144	0.704	191	181.693	136.31	1.333
6	5.794	9.437	0.614	162	5.853	6.257	0.935	192	79.848	63.049	1.266
7	8.577	12.718	0.674	163	6.43	4.973	1.293	193	240.989	145.212	1.660
8	4.491	4.176	1.075	164	0.395	1.705	0.232	194	185.607	176.692	1.050
9	12.77	17.11	0.746	165	2.532	2.717	0.932	195	126.013	110.683	1.139
10	8.794	11.648	0.755	166	4.739	5.78	0.820	196	102.756	118.139	0.870
11	10.107	14.921	0.677	167	4.141	6.733	0.615	197	167.541	180.644	0.927
12	7.155	9.388	0.762	168	5.758	7.456	0.772	198	71.611	56.052	1.278
13	8.329	13.555	0.614	169	0.316	0.319	0.991	199	521.659	456.524	1.143
14	7.353	9.699	0.758	170	0.383	0.858	0.446	200	327.406	439.069	0.746
15	10.103	13.573	0.744	171	7.887	12.687	0.622	201	123.286	99.666	1.237
		Mean	0.750			Mean	0.777			Mean	1.150
		SD	0.172			SD	0.281			SD	0.219
		%CV	22.923			%CV	36.105			%CV	19.014

CONCLUSIONS

Optimized for 10 µL Mitra® sampler procedure met the validation criteria with the use of Chromsystems kit reagents and materials.

Additional steps required for Mitra® sample preparation are simple, robust and automation friendly.

The VB/M correlations revealed that the capillary blood concentration of TAC and SIR are lower than in venous blood contrary to concentration of CSA. This data may support the fact that TAC and SIR enters erythrocytes in about 90%, while CSA only in 50% and the capillary blood has higher amount of serum fraction.

However, due to relatively small number of patients included in the study it needs further investigation on a larger group of patients and deeper statistical analysis